(90 × 210mm Size)

FUJIFILM

Wako

〈For Research Use Only〉 Code No. 294-65801 (1,000 Tests)

LabAssay™ Cholesterol (Cholesterol Oxidase·DAOS method)

[Introduction]

Cholesterol concentration in serum is known to be closely related with production, absorption and dissimilation in liver and intestinal tract, and lipoprotein metabolism in blood.

LabAssay™ Cholesterol is able to measure total cholesterol in mouse and human serum by an enzyme reaction using *N*-Ethyl-*N*-(2-hydroxy-3-sul-fopropyl)-3,5-dimethoxyaniline sodium salt (DAOS). It is a simultaneous multi-sample assay format using a microplate, but it measurements can also be made using a test tube.

[Kit contents]

(1)	Buffer (50mmol/L MES Buffer)	150mL	2vials
(2)	Chromogen After reconstitution, Cholesterol esterase 1.6units/mL Cholesterol oxidase 0.31units/mL Peroxidase (HRP) 5.2units/mL DAOS 0.95mmol/L 4-Aminoantipyrin 0.19mmol/L Ascorbate oxidase 4.4units/mL	for 150mL	2vials
(3)	Standard solution (Cholesterol 200mg/dL)	10mL	1vial

[Materials and apparatuses to be prepared]

- 96wells microplate (transparent type)
- Micropipette
- Incubator maintaining at 37°C*
- Plate mixer*
- Microplate reader with 600nm wavelength filter
(* if the microplate reader is not equipped.)

[for test tube method]

- Test tube
- Pipette
- Spectrophotometer or colorimeter with 600nm wavelength filter

[Assay principle]

Cholesterol esters in the sample are decomposed into free cholesterol and fatty acid by cholesterol esterase, when the sample reacts with chromo-gen reagent. The cholesterol is oxidized with existing free cholesterol by cholesterol oxidase, and simultaneously hydrogen peroxide is produced. The produced hydrogen peroxide let DAOS and 4-Aminoantipyrin oxidize and condensate quantitatively by peroxidase (HRP), which produces a blue pigment. Quantitation of total cholesterol in the sample can be made by measurement of the absorbance.

—1/8—

[Preparation of reagents to be used]

(1) Chromogen reagent :

Chromogen reagent is prepared by dissolving 1 vial of Chromogen with 1 vial of Buffer.

(2) Standard solution

Standard solution is prepared by dilution of the provided Standard.

No.	Provided Std.	Distilled water for dilution	Sample volume	Cholesterol concentration
1	50 μL	150 μL	2 μL	50.0mg/dL
2	100 μL	100 μL	2 μL	100.0mg/dL
3	Undiluted	—	2 μL	200.0mg/dL
4	Undiluted	—	4 μL	397.4mg/dL*
5	Undiluted	—	6 μL	592.2mg/dL*

*The test sample volume is usually 2 μL, but 4 or 6 μL are taken in this case. The cholesterol concentration must be corrected accordingly as indicated in the table above because of volume increasing.

[Procedure]

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 2 μL	Standard solution 2 μL	—
Chromogen reagent	300 μL	300 μL	300 μL
	Mix well, and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Main wavelength 600nm / Sub wavelength 700nm		


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 20 μL	Standard solution 20 μL	—
Chromogen reagent	3mL	3mL	3mL
	Mix well, and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600nm wavelength filter Spectrophotometer : Main wavelength 600nm (sub wavelength 700nm)		

—2/8—

[Standard curve]

Microplate reader : SAFIRE (TECAN)

[Performance]

- (1) Sensitivity
 - The absorbance is below 0.11, when measuring purified water as a sample.
 - The absorbance is 0.13 to 0.65, when measuring 200mg/dL cholesterol as a sample.
- (2) Specificity
 - The cholesterol concentration is less than $\pm 10\%$, when measuring the known concentration of control serum as a sample.

[Usage Notes]

- (1) Sample
 - Serum sample should be used immediately after collecting blood.
 - Hemolysis causes a slightly increase of the cholesterol value.
 - Ascorbic acid and bilirubin may not significantly affect the assay.
- (2) Notes on the assay
 - Measure the total cholesterol according to this insert sheet.
 - Use measurement apparatuses according to the use manual.
 - Keep the provided reagents under the indicated conditions before the expiration. Avoid using frozen reagents. Each of the reagents should be used immediately after opening, or store them under the indicated conditions.
 - Vials should be opened carefully, as the inner is decompressed.
 - The chromogen reaction is completed in ca. 2 min. Incubation up to 20 min. may not affect the assay. The reaction is stable and may not significantly affect the absorbance after 2 hours of the reaction.
 - This kit is for research use only, but not for diagnostic use.
- (3) Safety precautions
 - If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water.
 - A pipette should be used with a safety pipette filler.
 - Vials should be opened carefully.
- (4) Waste
 - The waste should be processed appropriately according to the law.

Expiration date : 1 year after the manufacture

Storage : Store at 2 ~ 10°C

Package : 1,000 tests

[Reference]

1. Allain, C. C., Poon, L. S., Chan, C. S. G., Richmond, W. and Fu, P. C. : *Clin. Chem.*, **20**, 470 (1974).
2. Richmond, W. : *Clin. Chem.*, **19**, 1350 (1973).
3. Richmond, W. : *Scand. J. Clin. Lab. Invest.*, **29** (Suppl.), 126 (1972).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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